

Peroxide test

0



1 mg/L H₂O₂



3



10



30



100

Instructions for use:

Dip test stick about 1 sec. into
test solution, shake off excess
liquid and read off at 5 sec.

Expiry date:
1 year after receipt

Keep cool
Contains 100 test strips



Roam Technology

Roam Technology NV

I.Z. Poort Genk 6835 - Geleenlaan 24 - 3600 Genk - Belgium

T +32 (0)89 44 00 42 - info@roamtechnology.com

www.roamtechnology.com

60620629